Veneuze beenulcera

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Referenties: 1. Lauchli, S., Bayard, I., Hafner, J. et al. [Healing times and the need for hospitalization for leg ulcers of different etiologies]. *Hautarzt* 2013; 64: 12, 917-922. 2. Hafner, J., Nobbe, S., Partsch, H. et al. Martorell hypertensive ischemic leg ulcer: a model of ischemic subcutaneous arteriosclerosis. *Arch Dermatol* 2010; 146: 9, 961-968. 3. Australian and New Zealand Clinical practice guidelines for prevention and management of Venous leg ulcers. 2011. 4. Hofman, D., Ryan, T.J., Arnold, F., et al. Pain in venous leg ulcers. *Journal of Wound Care* 1997; 6: 5, 222-224. [1997]. 5. Franks, P., Barker, J., Collier, M. et al. Management of patients with venous leg ulcer: challenges and current best practice. *J Wound Care*, 25; 6, Suppl, 1-67, 2016. 6. Wounds UK. Best Practice Statement: Holistic management of venous leg ulceration. 2016 London:Wounds UK. 7. Posnett, J., Gottrup, F., Lundgren, H., Saal, G. The resource impact of wounds on health-care providers in Europe. *J Wound Care* 2009;18: 4,154-61. 8. Abbade, L.P., Lastoria, S., de Almeida Rollo, H., et al. A sociodemographic, clinical study of patients with venous ulcer. *Int J Dermatol* 2005; 44: 12, 989-992. 9. McDaniel, H.B., Marston, W.A., Farber, M.A., et al. Recurrence of chronic venous ulcers on the basis of clinical, etiologic, anatomic, and pathophysiologic criteria and air plethysmography. *J Vasc Surg* 2002; 35: 4, 723-728. 10. Finlayson, K., Wu, M.L., Edwards, H.E. Identifying risk factors and protective factors for venous leg ulcer recurrence using a theoretical approach: A longitudinal study. *Int J Nurs Stud* 2015; 52: 6,1042-1045. 11. Harding, K., et al. Simplifying venous leg ulcer management. Consensus recommendations. *Wounds International* 2015. 12. Sufian, S., Lakhanpal, S., Marquez, J., et al. Superficial vein ablation for the treatment of primary chronic venous ulcers. *Phlebology* 2011; 26: 301-6. 13. O'Donnell, T.F., Passman, M.A., Marston, W.A. et al. Management of venous leg ulcers: Clinical practice guidelines of the Society for Vascular Surgery® and the American Venous Forum. *Journal of Vascular Surgery* 2014; 60: 2, 3S-59S. 14. Björk, R., Ehmman, S., S.T.R.I.D.E. Professional guide to compression garment selection for the lower extremity. *Journal of Wound Care* 2019; 28:16 suppl 1:1-44. 15. World Union of Wound Healing Societies (WUWHS). Principles of best practice: Compression in venous leg ulcers. A consensus document. London: MEP Ltd, 2008. 16. Humphreys, M., Stewart, A., Gohel, M. et al. Management of mixed arterial and venous leg ulcers. *Br J Surg* 2007; 94: 9, 1104-1147. 17. O'Meara, S., Cullum, N., Nelson, E.A., et al. Compression for venous leg ulcers. *Cochrane Database Syst Rev* 2012; 11: CD000265. 18. Nelson, E.A., Bell-Syer, S.E. Compression for preventing recurrence of venous ulcers. *Cochrane Database Syst Rev* 2012; 15: CD002303. 19. Wounds UK. Best Practice Statement: Ankle brachial pressure index (ABPI) in practice. London: Wounds UK, 2019. 20. Principles of compression in venous disease: a practitioner's guide to treatment and prevention of venous leg ulcers. *Wounds International*, 2013. 21. Wounds UK. Best Practice Statement: Addressing complexities in the management of venous leg ulcers. London: Wounds UK, 2019. 22. Gohel, M.S., Heatley, F., Liu, X. et al. A randomized trial of early endovenous ablation in venous ulceration. *N Engl J Med* 2018;378(22): 2105-14. 23. International Wound Infection Institute (IWII) Wound infection in clinical practice. *Wounds International* 2016. 24. Cutting, K.F., White, R. Defined and refined: criteria for identifying wound infection revisited. *British journal of community nursing* 2004;Mar;9(3): 6-15. 25. Gardener, S.E., Frantz, R.A., Doebbeling, B.N. The validity of the clinical signs and symptoms used to identify chronic wound infection, Wound repair and regeneration 2001;9: 178-186. 26. Bjarnsholt, T., Cooper, R., Fletcher, J. et al. World Union of Wound Healing Societies (WUWHS), Florence Congress, Position Document. Management of Biofilm. *Wounds International* 2016. 27. Swanson, T., Grothier, L., Schultz, G. Wound Infection Made Easy. *Wounds International* 2014. 28. Strohal, R., Gerber, V., Kröger, K., et al. Expert consensus on practical aspects of wound therapy with hemoglobin spray. *Wound management* 2016; 5: 276-284. 29. Teagle, A., Hargest, R. Management of pyoderma gangrenosum. *Journal of the Royal Society of Medicine* 2014;107(6):228-236. 30. White, R. Evidence for atraumatic soft silicone wound dressing use. *Wounds UK* 2005;1(3):104-109. 31. Waring, M., Biefeldt, S., Matzold, K.P., Butcher, M. An evaluation of the skin stripping of wound dressing adhesives. *J Wound Care* 2011;20:412-22. 32. White, R.A. Multinational survey of the assessment of pain when removing dressings. *Wounds UK* 2008;4:14-22. 33. David, F., Wutze, J.-L., Breton, N., et al. A randomised, controlled, non-inferiority trial comparing the performance of a soft silicone-coated wound contact layer (Mepitel One) with a lipidocolloid wound contact layer (UrgoTul) in the treatment of acute wounds. *Int Wound J* 2017 doi:10.1111/iwj.12853. 34. Gotschall, C.S., Morrison, M.I., Eichelberger, M.R. Prospective, randomized study of the efficacy of Mepitel on children with partial-thickness scalds. *J Burn Care Rehabil* 1998;19:279-83. 35. Silverstein, P., Heimbach, D., Meites, H., et al. An open, parallel, randomized, comparative, multicenter study to evaluate the cost effectiveness, performance, tolerance, and safety of a silver-containing soft silicone foam dressing (intervention) vs silver sulfadiazine cream. *J Burn Care Res* 2011;32:617-26. 36. Gee Kee, E.L., Kimble, R.M., Cuttle, L., Khan, A., Stockton, K.A. Randomized controlled trial of three burn dressings for partial thickness burns in children. *Burns* 2015;41:946-55. 37. Meaume, S., Van De Looverbosch, D., Heyman, H., Romanelli, M., Ciangherotti, A., Charpin, S. A study to compare a new self-adherent soft silicone dressing with a self-adherent polymer dressing in stage II pressure ulcers. *Ostomy Wound Manage* 2003; 49 (9): 44-51. 38. Wiberg, A.-B., Felli, F., Daun, E.-K. Preventing maceration with a soft silicone dressing: in vitro evaluation. Poster presentation at the 3rd Congress of the World Union of Wound Healing Societies, Toronto, Canada, 2008. 39. Bredow, J., Hoffmann, K., Hellmich, M., Eysel, P., Zarghooni, K. Randomized clinical trial to evaluate performance of flexible self-adherent absorbent dressing coated with silicone layer after hip, knee or spinal surgery in comparison to standard wound dressing. Poster presentation at the 5th Congress of the World Union of Wound Healing Societies, Florence, Italy, 2016. 40. Legsmatter:Org. Knees to toes: What you need to know. 41. Mölnlycke Health Care data on file. 42. Lantin, A., Diegel, C., Scheske, J., Schmitt, C., Brönnner, A., Jodl, H. Mepilex XT in practice: results of a study in German specialist wound care centres. *Wounds International* 2015; 6(4):18-22. 43. Mölnlycke Health Care. Mepilex® Border Flex Product Manual - Conformability PD-528870. Data on file. 44. Mölnlycke Health Care. Mepilex® Border Flex - External data - FEM simulations Mx-B Flex vs Allevyn Life. Report no. PD-529747. Data on file. 45. Mölnlycke Health Care. Mepilex® Border Flex-Bacteria encapsulation. Report no. PD-537072. Data on file. 46. Mölnlycke Health Care. Mepilex® Border Flex - Waterproofness. Report no. PD-532095. Data on file. 47. Chadwick, P., Taherinejad, F., Hamberg, K., Waring, M. Clinical and scientific data on a silver-containing soft-silicone foam dressing: an overview. *J Wound Care* 2009;18:483-9. 48. External lab report: NAMSAs 09C 29253 01/09C 29253 02. 49. Kles C.L., Murrach, C.P., Smith, K., et al. Achieving and sustaining zero. Preventing surgical site infections after isolated coronary artery bypass with saphenous vein harvest through implementation of a staff-driven quality improvement process. *Dimensions Crit Care Nurs* 2015;34:265-72. 50. Smet, S., Beele, H., Saine, L., Suys, E., Henrickx, B. Open, noncomparative, multi-centre post market clinical follow-up investigation to evaluate performance and safety on pressure ulcers when using a gelling fibre dressing as intend-ed. Poster Presentation at European Pressure Ulcer Advisory Panel Conference, 2015, Ghent, Belgium. 51. Chadwick, P., McCauley, J. Open, non-comparative, multicentre post clinical study of the performance and safety of a gelling fibre wound dressing on diabetic foot ulcers. *J Wound Care* 2016;25:290-300. 52. Davies, P., McCarty, S. An in-use product evaluation of a gelling fibre dressing in wound management. E-poster presentation at Wounds UK Conference, 2017, Harrogate, United Kingdom. 53. Lev-Tov et al. An interim analysis of clinical investigation to evaluate exudate management and comfort of use of an antimicrobial gelling fibre dressing* in medium to highly exudative wounds. Poster presented at the Symposium of Advanced Wound Care, Fall meeting 2018, Las Vegas, NV, USA. 54. Mölnlycke Health Care Laboratory Report PD-520425 (unpublished). Retention. 55. Mölnlycke Health Care. Exufiber. Gesellschaft für Versorgungskonzepte in der Wundbehandlung (GWW) mbH, Stuttgart, Germany. Data on file (unpublished report, 2017). 56. Mölnlycke Health Care. Data on file. [2018]. 57. Mölnlycke Health Care. Data on file. [2020]. 58. Mölnlycke Health Care Laboratory Report PD-521248 (unpublished). Tensile strength. 59. Mölnlycke Health Care Laboratory Report PD-521232 (unpublished). Absorption under compression. 60. Mölnlycke Health Care Laboratory Report PD-522900 (unpublished). Absorption of thick exudate. 61. Mölnlycke Health Care Laboratory Report PD-521245 (unpublished). Absorption of blood. 62. Hamberg, K., et al. Antimicrobial effect of a new silver-containing gelling fibre dressing against common wound pathogens. Poster presented at the Symposium on Advanced Wound Care Spring meeting/ Wound Healing Society (WHS) Annual Meeting 2017, Apr 05-09, 2017, San Diego, CA, USA. 63. Gil, J., et al. 2017. Evaluation of a Gelling Fiber Dressing with Silver to Eliminate MRSA Biofilm Infections and Enhance the Healing. Poster presented at the Symposium on Advanced Wound Care Spring meeting/ Wound Healing Society (WHS) Annual Meeting 2017, Apr 05-09, 2017, San Diego, CA, USA. 64. Davis, S. C., Li, J., Gil, J. et al. Preclinical evaluation of a novel silver gelling fibre dressing on *Pseudomonas aeruginosa* in a porcine wound infection model. *Wound Rep Reg* 2019;27: 360-365. 65. External lab report; NAMSAs 11C_51788_01. 66. Grocott, P. Clinical investigation Mepilex® Transfer, Clinical Investigation of a silicone dressing in product development phase in the palliative management of patients with pressure sores and malignant wounds, study id MIN101 UK, London UK, 2000. 67. Mölnlycke Health Care. Data on file, 2018. 68. Wiegand, C., et al. A superabsorbent polymer-containing wound dressing efficiently sequesters MMPs and inhibits collagenase activity in vitro. *Journal of Materials Science Materials in Medicine*. 2013;24(10):2473-2478. 69. Mölnlycke Health Care. Data on file, 2014. 70. Panca, M. C.K., Guest, J.F. Clinical and cost-effectiveness of absorbent dressings in the treatment of highly exuding VLU's. *J Wound Care*. 2013;22:109-8. 71. Mölnlycke Health Care. Data on file, 2019. 72. Mölnlycke Health Care. Data on file, 2016. 73. Hunt, S., Elg, F. Clinical effectiveness of hemoglobin spray (Granulox®) as adjunctive therapy in the treatment of chronic diabetic foot ulcers. *Diabetic Foot & Ankle* 2016; 7: 33101. 74. Elg, F., Hunt, S. Hemoglobin spray as adjunct therapy in complex wounds: Meta-analysis versus standard care alone in pooled data by wound type across three retrospective cohort controlled evaluations. *SAGE Open Med*. 2018 Jun; 6:2050312118784313. 75. Hunt, S., Elg, F. The clinical effectiveness of hemoglobin spray as adjunctive therapy in the treatment of chronic wounds. *J Wound Care*. 2017 Sep; 26(9):558-568. 76. *In vitro* suspension test [EN13727, EN 13624, EN 13704, EN 14476 - phase 2] with Granudacyn® wound irrigation solution. 77. Consensus on Wound Antisepsis: Update 2018, *Skin Pharmacol Physiol* 2018;31:28-58, DOI: 10.1159/000481545. 78. Harding, K.G., Vanscheidt, W., Partsch, H. et al. Adaptive compression therapy for venous leg ulcers: A clinically effective, patient-centred approach. *International Wound Journal*. 2016;13(3):317-25. 79. Eytier, C., Gazeau, E., Beneteau, G. et al. Convenience and tolerance of the combination of a soft silicone foam dressing and a two-way stretch tubular bandage in the management of local wounds. An observational study conducted by 304 registered nurses (RNs) on 2,401 patients. Results of TeMpO study. *Journal des plaies et cicatrisations*. 2013;88(18):38-44.